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- 1、Zhenyu Zheng,Taotao Zhou,Rong Hu,et al.A specific aggregation-induced emission-conjugated polymer enables visual monitoring of osteogenic differentiation.Bioactive Materials.July 2020.10.1016/j.bioactmat.2020.06.020. (IF 8.724)
- 2、Mei Li,Tenglong Li,Juan Yin,et al.Evaluation of toxicological effects of bisphenol S with an in vitro human bone marrow mesenchymal stem cell: Implications for bone health.TOXICOLOGY,December 2022.10.1016/j.tox.2022.153408. (IF 4.571)
- 3、Rui Mu,Bo Chen,Bo Bi,et al.LIM Mineralization Protein-1 Enhances the Committed Differentiation of Dental Pulp Stem Cells through the ERK1/2 and p38 MAPK Pathways and BMP Signaling.International Journal of Medical Sciences.July 2022.10.7150/ijms.70411. (IF 3.642)
- 4、Zhang Jie,Ma Zhenrong,Yan Kang,et al.Matrix Gla Protein Promotes the Bone Formation by Up-Regulating Wnt/ β -Catenin Signaling Pathway.Frontiers in Endocrinology,December 2019.10.3389/fendo.2019.00891. (IF 3.634)
- 5、Yu Fan,Dong Yan,Yang Yan-wei,et al.Effect of an Experimental Direct Pulp-capping Material on the Properties and Osteogenic Differentiation of Human Dental Pulp Stem Cells.Scientific Reports,October 2016.10.1038/srep34713.
- 6、Yuan Cao,Feng Chen,Suyue Zhu-Dongsheng,et al.Staphylococcus aureus infection initiates hypoxia-mediated STIP1 homology and U-box containing protein 1 upregulation to trigger osteomyelitis.TOXICON.July 2024.10.1016/j.toxicon.2024.108049.
- 7、Mahmoud M. Bakr,Mohamed Shamel,Shereen N. Raafat,et al.Effect of pulp capping materials on odontogenic differentiation of human dental pulp stem cells: An in vitro study.Clinical and Experimental Dental Research,December 2023.10.1002/cre2.816. (IF 1.8)

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