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- 1、Wu Baosheng,Xu Wenjie,Wu Kunjin,et al.Single-cell analysis of the amphioxus hepatic caecum and vertebrate liver reveals genetic mechanisms of vertebrate liver evolution.Nature Ecology & Evolution.August 2024.10.1038/s41559-024-02510-9.(IF 13.9)
- 2、Shijie Guo,Pengyu Wang,Yu Sun,et al.Transformation of Natural Resin Resina Draconis to 3D Functionalized Fibrous Scaffolds for Efficient Chronic Wound Healing.Advanced Healthcare Materials.June 2024.10.1002/adhm.202401105.(IF 10)
- 3、Yongqi Zhang,Jia Man,Jianing Liu,et al.Construction of the Mussel-Inspired PDAM/Lysine/Heparin Composite Coating Combining Multiple Anticoagulant Strategies.ACS Applied Materials & Interfaces.May 2023.10.1021/acsami.3c05077.(IF 9.5)
- 4、Jiaxi Liu,Xiang Zhou,Yin Zhang,et al.Rapid hemostasis and excellent antibacterial cerium-containing mesoporous bioactive glass/chitosan composite sponge for hemostatic material.Materials Today Chemistry.January 2022.10.1016/j.mtchem.2021.100735.(IF 7.613)
- 5、Yanmei He,Shimei Ruan,Guozhu Liang,et al.A Nonbactericidal Anionic Antimicrobial Peptide Provides Prophylactic and Therapeutic Efficacies against Bacterial Infections in Mice by Immunomodulatory–Antithrombotic Duality.JOURNAL OF MEDICINAL CHEMISTRY.April 2024.10.1021/acs.jmedchem.4c00342.(IF 6.8)

- 6、 Yumeng Wang,Kaicheng Du,Quanyou Wang,et al.A multidimensional strategy for characterization,distinction,and quality control of two Clinopodium medicinal plants.JOURNAL OF ETHNOPHARMACOLOGY.March 2024.10.1016/j.jep.2024.118019.(IF 4.8)
- 7、 Jingwen Liu,Lihua Geng,Jing Wang,et al.Structural Analysis and Anticoagulant Activity of Fucosylated Glycosaminoglycan from Sea Cucumber Phyllophorus proteus.Foods.September 2024.10.3390/foods13182889.(IF 4.7)
- 8、 Xinxin Pang,Yage Zhang,Zining Peng,et al.Hirudin reduces nephropathy microangiopathy in STZ-induced diabetes rats by inhibiting endothelial cell migration and angiogenesis.LIFE SCIENCES.May 2020.10.1016/j.lfs.2020.117779.(IF 3.647)
- 9、 Mingyi Shen,Changwei Zheng,Zhijin Hao,et al.Undescribed compounds from Clinopodium chinense (Benth.) O. Kuntze and their coagulation activity studies.FITOTERAPIA.November 2023.10.1016/j.fitote.2023.105736.(IF 3.4)
- 10、 Jiayi Liu,Xiang Zhou,Yin Zhang,et al.Rapid hemostasis and high bioactivity cerium-containing mesoporous bioglass for hemostatic materials.JOURNAL OF BIOMEDICAL MATERIALS RESEARCH PART B-APPLIED BIOMATERIALS.December 2021.10.1002/jbm.b.34996.(IF 3.368)
- 11、 Yong Zhao,Jiazhen Li,Fan Leng,et al.Degradable porous carboxymethyl chitin hemostatic microspheres.JOURNAL OF BIOMATERIALS SCIENCE-POLYMER EDITION.May 2020.10.1080/09205063.2020.1760461.(IF 2.69)
- 12、 Zhang Xiaojing,Mu Hongfang,Zhong Yan,et al.Effect of High Altitude Environment on Pharmacokinetic and Pharmacodynamic of Warfarin in Rats.CURRENT DRUG METABOLISM.February 2024.10.2174/0113892002277930240201101256.(IF 2.1)
- 13、 Nan-Nan Zhi,Kai Zong,Xiao-Yi Jia,et al.Effect of high pressure processing on fibrinolytic activity of fruit bromelain in vivo.JOURNAL OF FOOD PROCESS ENGINEERING.June 2019.10.1111/jfpe.13146.(IF 1.448)
- 14、 Ya-Qing Du,Ju Tang,Zhen-Xiang Zhang,et al.Correlation of Interleukin-18 and High-Sensitivity C-Reactive Protein with Perioperative Deep Vein Thrombosis in Patients with Ankle Fracture.ANNALS OF VASCULAR SURGERY.August 2018.10.1016/j.avsg.2018.06.013.(IF 1.363)

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